

# IMPLANT CATALOGUE



**PERTHESE™**  
BREAST IMPLANTS



**MENTOR®**

Make life more beautiful™

[www.mentormedical.co.uk](http://www.mentormedical.co.uk)  
Customer Service 0800 973 082  
[contact@medgb.jnj.com](mailto:contact@medgb.jnj.com)

NOT FOR DISTRIBUTION IN USA  
APPROVED FOR EMEA REGION  
© Mentor Worldwide LLC 2012 0907007 Rev C

Based on the high level of manufacturing standards and quality testing, Perouse Plastique was the first French implant manufacturer given the authorization to sell silicone gel products into the local market after a ban period of several years. The company launched its trademark **PERTHESE™ Breast Implant** in 1990.

During the following years, Perouse Plastique was not only able to continuously strengthen its market position in France but also in many other EU countries as well as in Latin America. With the acquisition through Mentor in 2007, more surgeons worldwide now have the opportunity to access this unique range of products adapted to the needs of medical professionals AND their patients in both Aesthetic and Reconstructive procedures.

The **PERTHESE™ Breast Implant** is a sterile gel-filled implant made of a trilaminar silicone envelope. The envelope consists of an internal and external layer made with a medical grade silicone elastomer and an intermediate barrier layer designed to reduce gel bleed.

The **PERTHESE™ Implant** range is available in two lines of breast implants:

### **PERTHESE™ Breast Implants**

Round breast implants filled with a soft and cohesive silicone gel, available either with a micro-textured surface or with a smooth surface.

Product Profile	Micro-textured	Smooth
TX, Low Profile	X	
MX, Moderate Profile	X	
MP, Moderate Plus Profile	X	
540, High Profile	X	
550, Ultra High Profile	X	
700T, Low Profile		X
510T, High Profile		X

### **PERTHESE™ Esthea Breast Implants**

Anatomical breast implants filled with a cohesive silicone gel with a micro-textured surface, available in different base shapes.

Product Profile	Micro-textured	Smooth
ELP-EHP-ESHP	X	
AX	X	
EFP	X	
EMP-EMHP	X	

All **PERTHESE™ Silicone Gel-filled Breast Implants** are covered by **PERTHESE Quality Assurance Program warranty** or **PatientCare Promise limited warranty†**

† In case of eligibility, depending on the date of surgery, the implant will fall under one of the above warranties.

## CONTENTS

### **PERTHESE™ Breast Implants Cohesive Silicone Gel Breast Implants, Round Shapes. . . . . 1**

PERTHESE™ Breast Implants TX, Low Profile, Micro-Textured  
PERTHESE™ Breast Implants MX, Moderate Profile, Micro-Textured  
PERTHESE™ Breast Implants MP, Moderate Plus Profile, Micro-Textured  
PERTHESE™ Breast Implants 540, High Profile, Micro-Textured  
PERTHESE™ Breast Implants 550, Ultra High Profile, Micro-Textured  
PERTHESE™ Breast Implants 700T, Low Profile, Smooth  
PERTHESE™ Breast Implants 510T, High Profile, Smooth  
PERTHESE™ Breast Implants Summary Chart

### **PERTHESE™ Esthea Breast Implants Cohesive Silicone Gel Breast Implants, Anatomical Shapes. . . . . 7**

PERTHESE™ Esthea Breast Implants ELP - EHP - ESHP, Micro-Textured  
PERTHESE™ Esthea Breast Implants AX, Micro-Textured  
PERTHESE™ Esthea Breast Implants EFP, Micro-Textured  
PERTHESE™ Esthea Breast Implants EMP - EMHP, Micro-Textured  
MENTOR® CPX4™ Tissue Expanders / PERTHESE™ Esthea Breast Implants Sizing Guide – Two-Stage Breast Reconstruction

### **PERTHESE™ Sterile Silicone Breast Implant Sizers, Round Shapes . . . . . 11**

PERTHESE™ Breast Implant Sizers TX, Low Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers MX, Moderate Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers 540, High Profile, Micro-Textured

### **PERTHESE™ Sterile Inflatable Breast Implant Sizers, Round Shapes . . . . . 13**

PERTHESE™ Breast Implant Sizers TX, Low Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers MX, Moderate Profile, Micro-Textured  
PERTHESE™ Breast Implant Sizers 540, High Profile, Micro-Textured

### **PERTHESE™ Esthea Sterile Silicone Breast Implant Sizers, Anatomical Shapes. . . . . 15**

PERTHESE™ Breast Esthea Implant Sizers ELP - EHP - ESHP, Micro-Textured  
PERTHESE™ Breast Esthea Implant Sizers AX, Micro-Textured  
PERTHESE™ Breast Esthea Implant Sizers EMP, Micro-Textured

### **PERTHESE™ Esthea Sterile Inflatable Breast Implant Sizers, Anatomical Shapes. . . . . 17**

PERTHESE™ Breast Implant Sizers ELP - EHP - ESHP, Micro-Textured  
PERTHESE™ Breast Implant Sizers AX, Micro-Textured  
PERTHESE™ Breast Implant Sizers EMP, Micro-Textured

### **Important Safety Information. . . . . 19**

### **Contact Information . . . . . 20**

Cohesive Silicone Gel Breast Implants,  
Round Shapes

PERTHESE™ Breast Implants TX,  
Low Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
100 cc	9,1 cm	2,2 cm	<b>TX 100</b>
125 cc	9,7 cm	2,4 cm	<b>TX 125</b>
150 cc	10,2 cm	2,6 cm	<b>TX 150</b>
170 cc	10,8 cm	2,7 cm	<b>TX 170</b>
185 cc	11,1 cm	2,8 cm	<b>TX 185</b>
200 cc	11,4 cm	2,9 cm	<b>TX 200</b>
220 cc	11,8 cm	2,9 cm	<b>TX 220</b>
240 cc	12,2 cm	3,0 cm	<b>TX 240</b>
260 cc	12,3 cm	3,1 cm	<b>TX 260</b>
280 cc	12,5 cm	3,2 cm	<b>TX 280</b>
300 cc	12,8 cm	3,2 cm	<b>TX 300</b>
325 cc	13,1 cm	3,3 cm	<b>TX 325</b>
350 cc	13,5 cm	3,4 cm	<b>TX 350</b>

PERTHESE™ Breast Implants MX,  
Moderate Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
100 cc	8,7 cm	2,8 cm	<b>MX 100</b>
125 cc	9,2 cm	3,0 cm	<b>MX 125</b>
150 cc	9,9 cm	3,2 cm	<b>MX 150</b>
175 cc	10,4 cm	3,3 cm	<b>MX 175</b>
200 cc	11,0 cm	3,4 cm	<b>MX 200</b>
225 cc	11,5 cm	3,5 cm	<b>MX 225</b>
250 cc	12,0 cm	3,6 cm	<b>MX 250</b>
275 cc	12,2 cm	3,7 cm	<b>MX 275</b>
300 cc	12,5 cm	3,8 cm	<b>MX 300</b>
325 cc	12,8 cm	4,0 cm	<b>MX 325</b>
350 cc	13,0 cm	4,2 cm	<b>MX 350</b>
375 cc	13,2 cm	4,4 cm	<b>MX 375</b>
400 cc	13,5 cm	4,5 cm	<b>MX 400</b>
450 cc	13,9 cm	4,6 cm	<b>MX 450</b>
500 cc	14,2 cm	5,0 cm	<b>MX 500</b>
550 cc	14,5 cm	5,2 cm	<b>MX 550</b>
600 cc	14,9 cm	5,4 cm	<b>MX 600</b>
650 cc	15,2 cm	5,5 cm	<b>MX 650</b>
700 cc	15,6 cm	5,6 cm	<b>MX 700</b>
800 cc	16,1 cm	6,0 cm	<b>MX 800</b>

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Cohesive Silicone Gel Breast Implants,  
Round Shapes CONT'D

PERTHESE™ Breast Implants MP,  
Moderate Plus Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
115 cc	8,7 cm	3,1 cm	<b>MP 115</b>
135 cc	9,2 cm	3,2 cm	<b>MP 135</b>
165 cc	9,9 cm	3,4 cm	<b>MP 165</b>
190 cc	10,5 cm	3,6 cm	<b>MP 190</b>
225 cc	11,0 cm	3,9 cm	<b>MP 225</b>
260 cc	11,6 cm	4,1 cm	<b>MP 260</b>
305 cc	12,1 cm	4,4 cm	<b>MP 305</b>
320 cc	12,3 cm	4,5 cm	<b>MP 320</b>
345 cc	12,5 cm	4,6 cm	<b>MP 345</b>
380 cc	12,7 cm	4,9 cm	<b>MP 380</b>
410 cc	13,0 cm	5,0 cm	<b>MP 410</b>
440 cc	13,2 cm	5,3 cm	<b>MP 440</b>
475 cc	13,5 cm	5,4 cm	<b>MP 475</b>
510 cc	13,8 cm	5,5 cm	<b>MP 510</b>
580 cc	14,3 cm	5,6 cm	<b>MP 580</b>
635 cc	14,4 cm	5,9 cm	<b>MP 635</b>
685 cc	14,9 cm	6,1 cm	<b>MP 685</b>
735 cc	15,2 cm	6,3 cm	<b>MP 735</b>
800 cc	15,7 cm	6,4 cm	<b>MP 800</b>
910 cc	16,2 cm	6,8 cm	<b>MP 910</b>

PERTHESE™ Breast Implants 540,  
High Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
75 cc	8,1 cm	2,3 cm	<b>540-075</b>
100 cc	8,2 cm	3,0 cm	<b>540-100</b>
125 cc	8,8 cm	3,2 cm	<b>540-125</b>
150 cc	9,2 cm	3,4 cm	<b>540-150</b>
175 cc	9,9 cm	3,6 cm	<b>540-175</b>
200 cc	10,1 cm	3,8 cm	<b>540-200</b>
230 cc	10,4 cm	4,0 cm	<b>540-230</b>
260 cc	11,0 cm	4,3 cm	<b>540-260</b>
280 cc	11,1 cm	4,4 cm	<b>540-280</b>
300 cc	11,2 cm	4,5 cm	<b>540-300</b>
325 cc	11,7 cm	4,6 cm	<b>540-325</b>
350 cc	11,9 cm	4,8 cm	<b>540-350</b>
380 cc	12,0 cm	4,9 cm	<b>540-380</b>
400 cc	12,1 cm	4,9 cm	<b>540-400</b>
440 cc	12,7 cm	5,3 cm	<b>540-440</b>
500 cc	13,2 cm	5,4 cm	<b>540-500</b>
550 cc	13,25 cm	5,9 cm	<b>540-550</b>
600 cc	13,8 cm	6,0 cm	<b>540-600</b>
650 cc	14,2 cm	6,1 cm	<b>540-650</b>
850 cc	15,5 cm	6,7 cm	<b>540-850</b>

## PERTHESE™ BREAST IMPLANTS

### Cohesive Silicone Gel Breast Implants, Round Shapes CONT'D

#### PERTHESE™ Breast Implants 550, Ultra High Profile, Micro-Textured



Volume	Diameter	Projection	Catalog #
95 cc	7,5 cm	3,5 cm	550-095
125 cc	8,2 cm	3,8 cm	550-125
150 cc	8,8 cm	3,9 cm	550-150
160 cc	9,0 cm	4,0 cm	550-160
175 cc	9,2 cm	4,1 cm	550-175
195 cc	9,5 cm	4,2 cm	550-195
215 cc	9,7 cm	4,3 cm	550-215
235 cc	10,1 cm	4,5 cm	550-235
265 cc	10,4 cm	4,7 cm	550-265
285 cc	10,8 cm	4,7 cm	550-285
315 cc	11,1 cm	4,9 cm	550-315
330 cc	11,2 cm	5,1 cm	550-330
345 cc	11,3 cm	5,2 cm	550-345
390 cc	11,7 cm	5,4 cm	550-390
410 cc	12,0 cm	5,4 cm	550-410
425 cc	12,1 cm	5,5 cm	550-425
445 cc	12,2 cm	5,7 cm	550-445
485 cc	12,6 cm	6,0 cm	550-485
530 cc	12,9 cm	6,0 cm	550-530
580 cc	13,4 cm	6,1 cm	550-580
630 cc	13,6 cm	6,3 cm	550-630
705 cc	14,1 cm	6,7 cm	550-705

### Cohesive Silicone Gel Breast Implants, Round Shapes CONT'D

#### PERTHESE™ Breast Implants 700T, Low Profile, Smooth



Volume	Diameter	Projection	Catalog #
120 cc	9,7 cm	2,4 cm	700 T-120
140 cc	10,3 cm	2,4 cm	700 T-140
160 cc	10,8 cm	2,5 cm	700 T-160
180 cc	11,2 cm	2,6 cm	700 T-180
200 cc	11,2 cm	2,8 cm	700 T-200
220 cc	11,9 cm	2,8 cm	700 T-220
240 cc	12,1 cm	2,9 cm	700 T-240
260 cc	12,2 cm	3,1 cm	700 T-260
280 cc	12,6 cm	3,2 cm	700 T-280
300 cc	12,8 cm	3,2 cm	700 T-300
325 cc	13,1 cm	3,3 cm	700 T-325
350 cc	13,5 cm	3,3 cm	700 T-350

#### PERTHESE™ Breast Implants 510T, High Profile, Smooth



Volume	Diameter	Projection	Catalog #
80 cc	8,1 cm	2,5 cm	510 T-080
100 cc	9,1 cm	2,7 cm	510 T-100
120 cc	9,3 cm	2,9 cm	510 T-120
140 cc	9,6 cm	3,3 cm	510 T-140
160 cc	9,6 cm	3,4 cm	510 T-160
180 cc	10,0 cm	3,6 cm	510 T-180
200 cc	10,2 cm	3,7 cm	510 T-200
220 cc	10,7 cm	3,9 cm	510 T-220
240 cc	10,7 cm	4,2 cm	510 T-240
270 cc	11,0 cm	4,3 cm	510 T-270
300 cc	11,4 cm	4,4 cm	510 T-300
340 cc	11,7 cm	4,6 cm	510 T-340
380 cc	11,9 cm	4,8 cm	510 T-380
400 cc	12,2 cm	5,0 cm	510 T-400
450 cc	12,6 cm	5,3 cm	510 T-450
500 cc	13,3 cm	5,5 cm	510 T-500


Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above



Cohesive Silicone Gel Breast Implants,  
Anatomical Shapes

Cohesive Silicone Gel Breast Implants,  
Anatomical Shapes CONT'D

LOCATION OF THE HIGHEST PROJECTION POINT (A) AND INFERIOR AND SUPERIOR CURVE LENGTH MEASUREMENTS



PERTHESE™ Esthea Breast Implants  
ELP - EHP - ESHP, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
190 cc	10,0 cm	12,5 cm	3,5 cm	3,6 cm	4,8 cm	9,2 cm	<b>ELP 190</b>
225 cc	10,0 cm	12,5 cm	4,0 cm	3,6 cm	5,6 cm	9,5 cm	<b>EHP 225</b>
235 cc	10,0 cm	12,5 cm	4,6 cm	3,5 cm	5,7 cm	9,5 cm	<b>ESHP 235</b>
240 cc	11,1 cm	13,8 cm	3,6 cm	3,7 cm	5,0 cm	10,2 cm	<b>ELP 240</b>
275 cc	11,1 cm	13,8 cm	4,2 cm	3,7 cm	5,8 cm	10,6 cm	<b>EHP 275</b>
295 cc	11,1 cm	13,8 cm	4,8 cm	3,7 cm	6,1 cm	10,6 cm	<b>ESHP 295</b>
270 cc	11,7 cm	14,0 cm	3,9 cm	4,1 cm	5,5 cm	10,3 cm	<b>ELP 270</b>
325 cc	11,7 cm	14,0 cm	4,7 cm	4,1 cm	6,4 cm	10,6 cm	<b>EHP 325</b>
355 cc	11,7 cm	14,0 cm	5,3 cm	4,0 cm	6,8 cm	10,7 cm	<b>ESHP 355</b>
310 cc	12,3 cm	14,4 cm	4,2 cm	4,3 cm	6,0 cm	10,4 cm	<b>ELP 310</b>
375 cc	12,3 cm	14,4 cm	5,1 cm	4,3 cm	6,8 cm	10,8 cm	<b>EHP 375</b>
415 cc	12,3 cm	14,4 cm	5,6 cm	4,3 cm	7,2 cm	10,9 cm	<b>ESHP 415</b>
360 cc	12,6 cm	14,8 cm	4,5 cm	4,7 cm	6,4 cm	10,5 cm	<b>ELP 360</b>
425 cc	12,6 cm	14,8 cm	5,4 cm	4,7 cm	7,4 cm	10,9 cm	<b>EHP 425</b>
475 cc	12,6 cm	14,8 cm	5,9 cm	4,7 cm	7,9 cm	11,1 cm	<b>ESHP 475</b>

PERTHESE™ Esthea Breast Implants  
AX, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
190 cc	10,6 cm	10,8 cm	4,3 cm	3,7 cm	5,8 cm	7,8 cm	<b>AX 190</b>
210 cc	10,9 cm	11,2 cm	4,4 cm	3,8 cm	6,0 cm	8,2 cm	<b>AX 210</b>
240 cc	11,5 cm	12,2 cm	4,6 cm	3,9 cm	6,1 cm	8,8 cm	<b>AX 240</b>
260 cc	11,9 cm	12,5 cm	4,8 cm	4,0 cm	6,3 cm	9,1 cm	<b>AX 260</b>
285 cc	12,1 cm	12,8 cm	4,9 cm	4,1 cm	6,4 cm	9,4 cm	<b>AX 285</b>
300 cc	12,6 cm	13,2 cm	4,9 cm	4,2 cm	6,6 cm	9,6 cm	<b>AX 300</b>
340 cc	12,9 cm	13,7 cm	4,9 cm	4,4 cm	6,9 cm	10,0 cm	<b>AX 340</b>
380 cc	13,3 cm	14,1 cm	5,0 cm	4,6 cm	7,2 cm	10,2 cm	<b>AX 380</b>
400 cc	13,6 cm	14,2 cm	5,1 cm	4,7 cm	7,4 cm	10,3 cm	<b>AX 400</b>
450 cc	13,9 cm	14,7 cm	5,2 cm	4,8 cm	7,5 cm	10,8 cm	<b>AX 450</b>
480 cc	14,5 cm	15,2 cm	5,3 cm	4,9 cm	7,7 cm	11,2 cm	<b>AX 480</b>
500 cc	14,9 cm	15,6 cm	5,3 cm	5,0 cm	7,9 cm	11,5 cm	<b>AX 500</b>

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

PERTHESE™ Esthea Breast Implants  
EFP, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
95 cc	8,1 cm	7,5 cm	3,4 cm	2,7 cm	4,2 cm	5,4 cm	<b>EFP 095</b>
115 cc	8,6 cm	8,1 cm	3,6 cm	3,0 cm	4,6 cm	5,8 cm	<b>EFP 115</b>
140 cc	9,2 cm	8,5 cm	3,8 cm	3,2 cm	5,0 cm	6,1 cm	<b>EFP 140</b>
180 cc	10,2 cm	9,2 cm	4,1 cm	3,5 cm	5,4 cm	6,6 cm	<b>EFP 180</b>
210 cc	10,7 cm	9,7 cm	4,3 cm	3,7 cm	5,7 cm	6,9 cm	<b>EFP 210</b>
250 cc	11,2 cm	10,2 cm	4,5 cm	3,9 cm	6,1 cm	7,3 cm	<b>EFP 250</b>
285 cc	11,7 cm	10,7 cm	4,7 cm	4,1 cm	6,4 cm	7,7 cm	<b>EFP 285</b>
325 cc	12,2 cm	11,3 cm	4,9 cm	4,2 cm	6,6 cm	8,1 cm	<b>EFP 325</b>
370 cc	12,7 cm	11,8 cm	5,1 cm	4,4 cm	6,9 cm	8,4 cm	<b>EFP 370</b>
415 cc	13,1 cm	12,2 cm	5,3 cm	4,7 cm	7,3 cm	8,7 cm	<b>EFP 415</b>
470 cc	13,6 cm	12,6 cm	5,5 cm	4,9 cm	7,7 cm	9,0 cm	<b>EFP 470</b>
520 cc	14,0 cm	13,1 cm	5,7 cm	5,1 cm	8,0 cm	9,4 cm	<b>EFP 520</b>
575 cc	14,6 cm	13,5 cm	5,9 cm	5,3 cm	8,3 cm	9,6 cm	<b>EFP 575</b>
640 cc	15,1 cm	14,0 cm	6,1 cm	5,6 cm	8,8 cm	9,9 cm	<b>EFP 640</b>
710 cc	15,5 cm	14,2 cm	6,3 cm	5,8 cm	9,1 cm	10,1 cm	<b>EFP 710</b>

PERTHESE™ Esthea Breast Implants  
EMP - EMHP, Micro-Textured



Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
80 cc	8,1 cm	6,8 cm	3,2 cm	2,4 cm	3,7 cm	4,9 cm	<b>EMP 080</b>
90 cc	8,1 cm	6,8 cm	3,5 cm	2,4 cm	4,1 cm	5,2 cm	<b>EMHP 090</b>
100 cc	8,6 cm	7,3 cm	3,5 cm	2,7 cm	4,2 cm	5,3 cm	<b>EMP 100</b>
110 cc	8,6 cm	7,3 cm	3,8 cm	2,7 cm	4,6 cm	5,5 cm	<b>EMHP 110</b>
120 cc	9,2 cm	7,7 cm	3,6 cm	2,9 cm	4,5 cm	5,5 cm	<b>EMP 120</b>
130 cc	9,2 cm	7,7 cm	4,0 cm	2,9 cm	4,9 cm	5,8 cm	<b>EMHP 130</b>
145 cc	10,2 cm	8,5 cm	3,9 cm	3,2 cm	5,0 cm	6,1 cm	<b>EMP 145</b>
170 cc	10,2 cm	8,5 cm	4,3 cm	3,1 cm	5,3 cm	6,4 cm	<b>EMHP 170</b>
165 cc	10,7 cm	9,1 cm	4,0 cm	3,4 cm	5,3 cm	6,4 cm	<b>EMP 165</b>
195 cc	10,7 cm	9,1 cm	4,5 cm	3,3 cm	5,6 cm	6,8 cm	<b>EMHP 195</b>
200 cc	11,2 cm	9,6 cm	4,3 cm	3,6 cm	5,7 cm	6,8 cm	<b>EMP 200</b>
235 cc	11,2 cm	9,6 cm	4,7 cm	3,5 cm	5,9 cm	7,2 cm	<b>EMHP 235</b>
230 cc	11,7 cm	10,1 cm	4,5 cm	3,8 cm	6,0 cm	7,2 cm	<b>EMP 230</b>
270 cc	11,7 cm	10,1 cm	4,9 cm	3,7 cm	6,3 cm	7,6 cm	<b>EMHP 270</b>
265 cc	12,2 cm	10,6 cm	4,6 cm	3,9 cm	6,1 cm	7,6 cm	<b>EMP 265</b>
320 cc	12,2 cm	10,6 cm	5,1 cm	3,9 cm	6,6 cm	7,9 cm	<b>EMHP 320</b>
305 cc	12,7 cm	11,1 cm	4,9 cm	4,1 cm	6,4 cm	7,9 cm	<b>EMP 305</b>
360 cc	12,7 cm	11,1 cm	5,3 cm	4,1 cm	6,9 cm	8,3 cm	<b>EMHP 360</b>
345 cc	13,1 cm	11,6 cm	5,0 cm	4,4 cm	6,8 cm	8,2 cm	<b>EMP 345</b>
405 cc	13,1 cm	11,6 cm	5,5 cm	4,3 cm	7,2 cm	8,6 cm	<b>EMHP 405</b>
390 cc	13,6 cm	12,0 cm	5,2 cm	4,6 cm	7,2 cm	8,5 cm	<b>EMP 390</b>
455 cc	13,6 cm	12,0 cm	5,7 cm	4,5 cm	7,6 cm	9,0 cm	<b>EMHP 455</b>
435 cc	14,0 cm	12,5 cm	5,5 cm	4,8 cm	7,5 cm	8,9 cm	<b>EMP 435</b>
490 cc	14,0 cm	12,5 cm	5,8 cm	4,3 cm	7,7 cm	9,6 cm	<b>EMHP 490</b>
475 cc	14,6 cm	13,0 cm	5,7 cm	4,9 cm	7,7 cm	9,2 cm	<b>EMP 475</b>
545 cc	14,6 cm	13,0 cm	6,0 cm	4,6 cm	8,0 cm	9,7 cm	<b>EMHP 545</b>
535 cc	15,1 cm	13,3 cm	5,9 cm	5,0 cm	8,1 cm	9,6 cm	<b>EMP 535</b>
625 cc	15,1 cm	13,3 cm	6,2 cm	5,0 cm	8,6 cm	9,9 cm	<b>EMHP 625</b>
585 cc	15,5 cm	13,6 cm	6,1 cm	5,2 cm	8,4 cm	9,8 cm	<b>EMP 585</b>
680 cc	15,5 cm	13,6 cm	6,4 cm	5,2 cm	8,9 cm	10,2 cm	<b>EMHP 680</b>





## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Silicone Gel Breast Implant Sizers, Round Shapes

#### PERTHESE™ Breast Implants TX, Sterile Gel Sizers\* Low Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
TX 100	100 cc	9,1 cm	2,2 cm	GABTX100STE
TX 125	125 cc	9,7 cm	2,4 cm	GABTX125STE
TX 150	150 cc	10,2 cm	2,6 cm	GABTX150STE
TX 170	170 cc	10,8 cm	2,7 cm	GABTX170STE
TX 185	185 cc	11,1 cm	2,8 cm	GABTX185STE
TX 200	200 cc	11,4 cm	2,9 cm	GABTX200STE
TX 220	220 cc	11,8 cm	2,9 cm	GABTX220STE
TX 240	240 cc	12,2 cm	3,0 cm	GABTX240STE
TX 260	260 cc	12,3 cm	3,1 cm	GABTX260STE
TX 280	280 cc	12,5 cm	3,2 cm	GABTX280STE
TX 300	300 cc	12,8 cm	3,2 cm	GABTX300STE
TX 325	325 cc	13,1 cm	3,3 cm	GABTX325STE
TX 350	350 cc	13,5 cm	3,4 cm	GABTX350STE

#### PERTHESE™ Breast Implants MX, Sterile Gel Sizers\* Moderate Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
MX 100	100 cc	8,7 cm	2,8 cm	GABMX100STE
MX 125	125 cc	9,2 cm	3,0 cm	GABMX125STE
MX 150	150 cc	9,9 cm	3,2 cm	GABMX150STE
MX 175	175 cc	10,4 cm	3,3 cm	GABMX175STE
MX 200	200 cc	11,0 cm	3,4 cm	GABMX200STE
MX 225	225 cc	11,5 cm	3,5 cm	GABMX225STE
MX 250	250 cc	12,0 cm	3,6 cm	GABMX250STE
MX 275	275 cc	12,2 cm	3,7 cm	GABMX275STE
MX 300	300 cc	12,5 cm	3,8 cm	GABMX300STE
MX 325	325 cc	12,8 cm	4,0 cm	GABMX325STE
MX 350	350 cc	13,0 cm	4,2 cm	GABMX350STE
MX 375	375 cc	13,2 cm	4,4 cm	GABMX375STE
MX 400	400 cc	13,5 cm	4,5 cm	GABMX400STE
MX 450	450 cc	13,9 cm	4,6 cm	GABMX450STE
MX 500	500 cc	14,2 cm	5,0 cm	GABMX500STE
MX 550	550 cc	14,5 cm	5,2 cm	GABMX550STE
MX 600	600 cc	14,9 cm	5,4 cm	GABMX600STE
MX 650	650 cc	15,2 cm	5,5 cm	GABMX650STE
MX 700	700 cc	15,6 cm	5,6 cm	GABMX700STE
MX 800	800 cc	16,1 cm	6,0 cm	GABMX800STE

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Silicone Gel Breast Implant Sizers, Round Shapes CONT'D

#### PERTHESE™ Breast Implants 540, Sterile Gel Sizers\* High Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
540-075	75 cc	8,1 cm	2,3 cm	GAB540075STE
540-100	100 cc	8,2 cm	3,0 cm	GAB540100STE
540-125	125 cc	8,8 cm	3,2 cm	GAB540125STE
540-150	150 cc	9,2 cm	3,4 cm	GAB540150STE
540-175	175 cc	9,9 cm	3,6 cm	GAB540175STE
540-200	200 cc	10,1 cm	3,8 cm	GAB540200STE
540-230	230 cc	10,4 cm	4,0 cm	GAB540230STE
540-260	260 cc	11,0 cm	4,3 cm	GAB540260STE
540-280	280 cc	11,1 cm	4,4 cm	GAB540280STE
540-300	300 cc	11,2 cm	4,5 cm	GAB540300STE
540-325	325 cc	11,7 cm	4,6 cm	GAB540325STE
540-350	350 cc	11,9 cm	4,8 cm	GAB540350STE
540-380	380 cc	12,0 cm	4,9 cm	GAB540380STE
540-400	400 cc	12,1 cm	4,9 cm	GAB540400STE
540-440	440 cc	12,7 cm	5,3 cm	GAB540440STE
540-500	500 cc	13,2 cm	5,4 cm	GAB540500STE
540-550	550 cc	13,25 cm	5,9 cm	GAB540550STE
540-600	600 cc	13,8 cm	6,0 cm	GAB540600STE
540-650	650 cc	14,2 cm	6,1 cm	GAB540650STE
540-850	850 cc	15,5 cm	6,7 cm	GAB540850STE

\*Sterile - Single Use Only

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Inflatable Sizers, Round Shapes

#### PERTHESE™ Breast Implants TX, Sterile Inflatable Sizers\* Low Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
TX 100	100 cc	9,1 cm	2,2 cm	GGSTX100
TX 125	125 cc	9,7 cm	2,4 cm	GGSTX125
TX 150	150 cc	10,2 cm	2,6 cm	GGSTX150
TX 170	170 cc	10,8 cm	2,7 cm	GGSTX170
TX 185	185 cc	11,1 cm	2,8 cm	GGSTX185
TX 200	200 cc	11,4 cm	2,9 cm	GGSTX200
TX 220	220 cc	11,8 cm	2,9 cm	GGSTX220
TX 240	240 cc	12,2 cm	3,0 cm	GGSTX240
TX 260	260 cc	12,3 cm	3,1 cm	GGSTX260
TX 280	280 cc	12,5 cm	3,2 cm	GGSTX280
TX 300	300 cc	12,8 cm	3,2 cm	GGSTX300
TX 325	325 cc	13,1 cm	3,3 cm	GGSTX325
TX 350	350 cc	13,5 cm	3,4 cm	GGSTX350

#### PERTHESE™ Breast Implants MX, Sterile Inflatable Sizers\* Moderate Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
MX 125	125 cc	9,2 cm	3,0 cm	GGSMX125
MX 150	150 cc	9,9 cm	3,2 cm	GGSMX150
MX 175	175 cc	10,4 cm	3,3 cm	GGSMX175
MX 200	200 cc	11,0 cm	3,4 cm	GGSMX200
MX 225	225 cc	11,5 cm	3,5 cm	GGSMX225
MX 250	250 cc	12,0 cm	3,6 cm	GGSMX250
MX 275	275 cc	12,2 cm	3,7 cm	GGSMX275
MX 300	300 cc	12,5 cm	3,8 cm	GGSMX300
MX 325	325 cc	12,8 cm	4,0 cm	GGSMX325
MX 350	350 cc	13,0 cm	4,2 cm	GGSMX350
MX 375	375 cc	13,2 cm	4,4 cm	GGSMX375
MX 400	400 cc	13,5 cm	4,5 cm	GGSMX400

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

## PERTHESE™ BREAST IMPLANT SIZERS

### Sterile Inflatable Sizers, Round Shapes CONT'D

#### PERTHESE™ Breast Implants 540, Sterile Inflatable Sizers\* High Profile, Micro-Textured



Definitive Implant	Volume	Diameter	Projection	Catalog #
540-200	200 cc	10,1 cm	3,8 cm	GG540200
540-230	230 cc	10,4 cm	4,0 cm	GG540230
540-260	260 cc	11,0 cm	4,3 cm	GG540260
540-280	280 cc	11,1 cm	4,4 cm	GG540280
540-300	300 cc	11,2 cm	4,5 cm	GG540300
540-325	325 cc	11,7 cm	4,6 cm	GG540325
540-350	350 cc	11,9 cm	4,8 cm	GG540350
540-380	380 cc	12,0 cm	4,9 cm	GG540380
540-400	400 cc	12,1 cm	4,9 cm	GG540400

\*Sterile - Single Use Only

Sterile Silicone Gel Breast Implant Sizers,  
Anatomical Shapes

PERTHESE™ Esthea Breast Implants, Sterile Gel Sizers\*  
ELP - EHP - ESHP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
ELP 190	190 cc	10,0 cm	12,5 cm	3,5 cm	3,6 cm	4,8 cm	9,2 cm	GABELP190STE
EHP 225	225 cc	10,0 cm	12,5 cm	4,0 cm	3,6 cm	5,6 cm	9,5 cm	GABEHP225STE
ESHP 235	235 cc	10,0 cm	12,5 cm	4,6 cm	3,5 cm	5,7 cm	9,5 cm	GABESH235STE
ELP 240	240 cc	11,1 cm	13,8 cm	3,6 cm	3,7 cm	5,0 cm	10,2 cm	GABELP240STE
EHP 275	275 cc	11,1 cm	13,8 cm	4,2 cm	3,7 cm	5,8 cm	10,6 cm	GABEHP275STE
ESHP 295	295 cc	11,1 cm	13,8 cm	4,8 cm	3,7 cm	6,1 cm	10,6 cm	GABESH295STE
ELP 270	270 cc	11,7 cm	14,0 cm	3,9 cm	4,1 cm	5,5 cm	10,3 cm	GABELP270STE
EHP 325	325 cc	11,7 cm	14,0 cm	4,7 cm	4,1 cm	6,4 cm	10,6 cm	GABEHP325STE
ESHP 355	355 cc	11,7 cm	14,0 cm	5,3 cm	4,0 cm	6,8 cm	10,7 cm	GABESH355STE
ELP 310	310 cc	12,3 cm	14,4 cm	4,2 cm	4,3 cm	6,0 cm	10,4 cm	GABELP310STE
EHP 375	375 cc	12,3 cm	14,4 cm	5,1 cm	4,3 cm	6,8 cm	10,8 cm	GABEHP375STE
ESHP 415	415 cc	12,3 cm	14,4 cm	5,6 cm	4,3 cm	7,2 cm	10,9 cm	GABESH415STE
ELP 360	360 cc	12,6 cm	14,8 cm	4,5 cm	4,7 cm	6,4 cm	10,5 cm	GABELP360STE
EHP 425	425 cc	12,6 cm	14,8 cm	5,4 cm	4,7 cm	7,4 cm	10,9 cm	GABEHP425STE
ESHP 475	475 cc	12,6 cm	14,8 cm	5,9 cm	4,7 cm	7,9 cm	11,1 cm	GABESH475STE

PERTHESE™ Esthea Breast Implants, Sterile Gel Sizers\*  
AX, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
AX 190	190 cc	10,6 cm	10,8 cm	4,3 cm	3,7 cm	5,8 cm	7,8 cm	GABAX190STE
AX 210	210 cc	10,9 cm	11,2 cm	4,4 cm	3,8 cm	6,0 cm	8,2 cm	GABAX210STE
AX 240	240 cc	11,5 cm	12,2 cm	4,6 cm	3,9 cm	6,1 cm	8,8 cm	GABAX240STE
AX 260	260 cc	11,9 cm	12,5 cm	4,8 cm	4,0 cm	6,3 cm	9,1 cm	GABAX260STE
AX 285	285 cc	12,1 cm	12,8 cm	4,9 cm	4,1 cm	6,4 cm	9,4 cm	GABAX285STE
AX 300	300 cc	12,6 cm	13,2 cm	4,9 cm	4,2 cm	6,6 cm	9,6 cm	GABAX300STE
AX 340	340 cc	12,9 cm	13,7 cm	4,9 cm	4,4 cm	6,9 cm	10,0 cm	GABAX340STE
AX 380	380 cc	13,3 cm	14,1 cm	5,0 cm	4,6 cm	7,2 cm	10,2 cm	GABAX380STE
AX 400	400 cc	13,6 cm	14,2 cm	5,1 cm	4,7 cm	7,4 cm	10,3 cm	GABAX400STE
AX 450	450 cc	13,9 cm	14,7 cm	5,2 cm	4,8 cm	7,5 cm	10,8 cm	GABAX450STE
AX 480	480 cc	14,5 cm	15,2 cm	5,3 cm	4,9 cm	7,7 cm	11,2 cm	GABAX480STE
AX 500	500 cc	14,9 cm	15,6 cm	5,3 cm	5,0 cm	7,9 cm	11,5 cm	GABAX500STE

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.  
Not all units will conform exactly to the dimensions noted above

Sterile Silicone Gel Breast Implant Sizers,  
Anatomical Shapes CONT'D

PERTHESE™ Esthea Breast Implants, Sterile Gel Sizers\*  
EMP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
EMP 080	80 cc	8,1 cm	6,8 cm	3,2 cm	2,4 cm	3,7 cm	4,9 cm	GABEMP080STE
EMP 100	100 cc	8,6 cm	7,3 cm	3,5 cm	2,7 cm	4,2 cm	5,3 cm	GABEMP100STE
EMP 120	120 cc	9,2 cm	7,7 cm	3,6 cm	2,9 cm	4,5 cm	5,5 cm	GABEMP120STE
EMP 145	145 cc	10,2 cm	8,5 cm	3,9 cm	3,2 cm	5,0 cm	6,1 cm	GABEMP145STE
EMP 165	165 cc	10,7 cm	9,1 cm	4,0 cm	3,4 cm	5,3 cm	6,4 cm	GABEMP165STE
EMP 200	200 cc	11,2 cm	9,6 cm	4,3 cm	3,6 cm	5,7 cm	6,8 cm	GABEMP200STE
EMP 230	230 cc	11,7 cm	10,1 cm	4,5 cm	3,8 cm	6,0 cm	7,2 cm	GABEMP230STE
EMP 265	265 cc	12,2 cm	10,6 cm	4,6 cm	3,9 cm	6,1 cm	7,6 cm	GABEMP265STE
EMP 305	305 cc	12,7 cm	11,1 cm	4,9 cm	4,1 cm	6,4 cm	7,9 cm	GABEMP305STE
EMP 345	345 cc	13,1 cm	11,6 cm	5,0 cm	4,4 cm	6,8 cm	8,2 cm	GABEMP345STE
EMP 390	390 cc	13,6 cm	12,0 cm	5,2 cm	4,6 cm	7,2 cm	8,5 cm	GABEMP390STE
EMP 435	435 cc	14,0 cm	12,5 cm	5,5 cm	4,8 cm	7,5 cm	8,9 cm	GABEMP435STE
EMP 475	475 cc	14,6 cm	13,0 cm	5,7 cm	4,9 cm	7,7 cm	9,2 cm	GABEMP475STE
EMP 535	535 cc	15,1 cm	13,3 cm	5,9 cm	5,0 cm	8,1 cm	9,6 cm	GABEMP535STE
EMP 585	585 cc	15,5 cm	13,6 cm	6,1 cm	5,2 cm	8,4 cm	9,8 cm	GABEMP585STE

\*Sterile - Single Use Only

Sterile Inflatable Sizers,  
Anatomical Shapes

PERTHESE™ Esthea Breast Implants, Sterile Inflatable Sizers\*  
ELP - EHP - ESHP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
ELP 190	190 cc	10,0 cm	12,5 cm	3,5 cm	3,6 cm	4,8 cm	9,2 cm	GGSELP190
EHP 225	225 cc	10,0 cm	12,5 cm	4,0 cm	3,6 cm	5,6 cm	9,5 cm	GGSEHP225
ESHP 235	235 cc	10,0 cm	12,5 cm	4,6 cm	3,5 cm	5,7 cm	9,5 cm	GGSESH235
ELP 240	240 cc	11,1 cm	13,8 cm	3,6 cm	3,7 cm	5,0 cm	10,2 cm	GGSELP240
EHP 275	275 cc	11,1 cm	13,8 cm	4,2 cm	3,7 cm	5,8 cm	10,6 cm	GGSEHP275
ESHP 295	295 cc	11,1 cm	13,8 cm	4,8 cm	3,7 cm	6,1 cm	10,6 cm	GGSESH295
ELP 270	270 cc	11,7 cm	14,0 cm	3,9 cm	4,1 cm	5,5 cm	10,3 cm	GGSELP270
EHP 325	325 cc	11,7 cm	14,0 cm	4,7 cm	4,1 cm	6,4 cm	10,6 cm	GGSEHP325
ESHP 355	355 cc	11,7 cm	14,0 cm	5,3 cm	4,0 cm	6,8 cm	10,7 cm	GGSESH355
ELP 310	310 cc	12,3 cm	14,4 cm	4,2 cm	4,3 cm	6,0 cm	10,4 cm	GGSELP310
EHP 375	375 cc	12,3 cm	14,4 cm	5,1 cm	4,3 cm	6,8 cm	10,8 cm	GGSEHP375
ESHP 415	415 cc	12,3 cm	14,4 cm	5,6 cm	4,3 cm	7,2 cm	10,9 cm	GGSESH415
ELP 360	360 cc	12,6 cm	14,8 cm	4,5 cm	4,7 cm	6,4 cm	10,5 cm	GGSELP360
EHP 425	425 cc	12,6 cm	14,8 cm	5,4 cm	4,7 cm	7,4 cm	10,9 cm	GGSEHP425
ESHP 475	475 cc	12,6 cm	14,8 cm	5,9 cm	4,7 cm	7,9 cm	11,1 cm	GGSESH475

PERTHESE™ Esthea Breast Implants, Sterile Inflatable Sizers\*  
AX, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
AX 240	240 cc	11,5 cm	12,2 cm	4,6 cm	3,9 cm	6,1 cm	8,8 cm	GGSA240
AX 260	260 cc	11,9 cm	12,5 cm	4,8 cm	4,0 cm	6,3 cm	9,1 cm	GGSA260
AX 285	285 cc	12,1 cm	12,8 cm	4,9 cm	4,1 cm	6,4 cm	9,4 cm	GGSA285
AX 300	300 cc	12,6 cm	13,2 cm	4,9 cm	4,2 cm	6,6 cm	9,6 cm	GGSA300
AX 340	340 cc	12,9 cm	13,7 cm	4,9 cm	4,4 cm	6,9 cm	10,0 cm	GGSA340
AX 380	380 cc	13,3 cm	14,1 cm	5,0 cm	4,6 cm	7,2 cm	10,2 cm	GGSA380
AX 400	400 cc	13,6 cm	14,2 cm	5,1 cm	4,7 cm	7,4 cm	10,3 cm	GGSA400

\*Sterile - Single Use Only

Note: Individual implant dimensions may vary slightly in products of this type.

Not all units will conform exactly to the dimensions noted above

Sterile Inflatable Sizers,  
Anatomical Shapes CONT'D

PERTHESE™ Esthea Breast Implants, Sterile Inflatable Sizers\*  
EMP, Micro-Textured



Definitive Implant	Volume	Width	Height	Projection	A	Inferior Curve Length	Superior Curve Length	Catalog #
EMP 120	120 cc	9,2 cm	7,7 cm	3,6 cm	2,9 cm	4,5 cm	5,5 cm	GGSEMP120
EMP 145	145 cc	10,2 cm	8,5 cm	3,9 cm	3,2 cm	5,0 cm	6,1 cm	GGSEMP145
EMP 165	165 cc	10,7 cm	9,1 cm	4,0 cm	3,4 cm	5,3 cm	6,4 cm	GGSEMP165
EMP 200	200 cc	11,2 cm	9,6 cm	4,3 cm	3,6 cm	5,7 cm	6,8 cm	GGSEMP200
EMP 230	230 cc	11,7 cm	10,1 cm	4,5 cm	3,8 cm	6,0 cm	7,2 cm	GGSEMP230
EMP 265	265 cc	12,2 cm	10,6 cm	4,6 cm	3,9 cm	6,1 cm	7,6 cm	GGSEMP265
EMP 305	305 cc	12,7 cm	11,1 cm	4,9 cm	4,1 cm	6,4 cm	7,9 cm	GGSEMP305
EMP 345	345 cc	13,1 cm	11,6 cm	5,0 cm	4,4 cm	6,8 cm	8,2 cm	GGSEMP345
EMP 390	390 cc	13,6 cm	12,0 cm	5,2 cm	4,6 cm	7,2 cm	8,5 cm	GGSEMP390
EMP 435	435 cc	14,0 cm	12,5 cm	5,5 cm	4,8 cm	7,5 cm	8,9 cm	GGSEMP435
EMP 475	475 cc	14,6 cm	13,0 cm	5,7 cm	4,9 cm	7,7 cm	9,2 cm	GGSEMP475
EMP 535	535 cc	15,1 cm	13,3 cm	5,9 cm	5,0 cm	8,1 cm	9,6 cm	GGSEMP535
EMP 585	585 cc	15,5 cm	13,6 cm	6,1 cm	5,2 cm	8,4 cm	9,8 cm	GGSEMP585

\*Sterile - Single Use Only

## Important Safety Information

MENTOR® Breast Implants are indicated for breast augmentation, in women who are at least 18 years old, or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body, with existing cancer or pre-cancer of their breast(s) who have not received adequate treatment for those conditions or who are pregnant or nursing.

There are risks associated with breast implant surgery. Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery.

The most common complications with MENTOR® MemoryGel™ Breast Implants and PERTHESE™ Breast Implants include re-operation, implant removal, capsular contracture, asymmetry, and breast pain. A lower risk of complication is implant rupture, which is most often silent. The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture.

The most common complications with MENTOR® Saline-Filled Implants include re-operation, implant removal, capsular contracture, wrinkling, deflation, asymmetry, and breast pain.

MENTOR® CONTOUR PROFILE™ Tissue Expanders and MENTOR® Smooth Round Tissue Expanders with Remote Dome are used for breast reconstruction following mastectomy. These expanders are intended for temporary subcutaneous or submuscular implantation and are not intended for use beyond six months. CONTOUR PROFILE™ Tissue Expanders are devices that contain magnetic injection domes and are NOT MRI compatible. Do not use the CONTOUR PROFILE™ Tissue Expander in patients where an MRI may be needed. DO NOT use the CONTOUR PROFILE™ Tissue Expander in patients that have a previously implanted device that could be affected by a magnetic field. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed.

Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery.

For detailed indications, contraindications, warning and precautions associated with the use of all MENTOR® Implantable Devices, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at [www.mentorwllc.eu](http://www.mentorwllc.eu)

## Our Credo

We believe our first responsibility is to the doctors, nurses and patients, to mothers and fathers and all others who use our products and services. In meeting their needs everything we do must be of high quality. We must constantly strive to reduce our costs in order to maintain reasonable prices. Customers' orders must be serviced promptly and accurately. Our suppliers and distributors must have an opportunity to make a fair profit.

We are responsible to our employees, the men and women who work with us throughout the world. Everyone must be considered as an individual. We must respect their dignity and recognize their merit. They must have a sense of security in their jobs. Compensation must be fair and adequate, and working conditions clean, orderly and safe. We must be mindful of ways to help our employees fulfill their family responsibilities. Employees must feel free to make suggestions and complaints. There must be equal opportunity for employment, development and advancement for those qualified. We must provide competent management, and their actions must be just and ethical.

We are responsible to the communities in which we live and work and to the world community as well. We must be good citizens – support good works and charities and bear our fair share of taxes. We must encourage civic improvements and better health and education. We must maintain in good order the property we are privileged to use, protecting the environment and natural resources.

Our final responsibility is to our stockholders. Business must make a sound profit. We must experiment with new ideas. Research must be carried on, innovative programs developed and mistakes paid for. New equipment must be purchased, new facilities provided and new products launched. Reserves must be created to provide for adverse times. When we operate according to these principles, the stockholders should realize a fair return.



A division of Johnson & Johnson Medical Limited  
Pinewood Campus, Nine Mile Ride, Wokingham  
Berkshire, RG40 3EW  
Customer Service 0800 973 082  
[contact@medgb.jnj.com](mailto:contact@medgb.jnj.com)  
[www.mentormedical.co.uk](http://www.mentormedical.co.uk)